


KATANA™ ZIRCONIA RESTORATION SYSTEM APPROPRIATE MATCHING MATERIALS



BENEFITS IN A DAILY WORKFLOW

- ✓ Developed and produced in Japan ensuring high product quality
- ✓ Natural color gradient within multi-layered discs as a foundation for true-to-life restorations
- ✓ Trueness to VITA classical A1-D4 shade guide
- ✓ Great mechanical properties allowing for less invasive preparations and more complex posterior restorations than glass ceramics
- ✓ Compatibility with self-adhesive cements supporting simplified cementation protocols




Translucency / Flexural Strength	Laminate veneer Inlay/Onlay	Anterior Crown	Posterior Crown	3-unit Bridge	Long-span Bridge	Crown/ Bridges for implants
KATANA™ Zirconia UTML 43% - 557 MPa	✓	✓	✓	✓*		✓
KATANA™ Zirconia STML 38% - 748 MPa		✓	✓	✓		✓
KATANA™ Zirconia HTML™ 31% - 1125 MPa		✓	✓	✓	✓	✓
KATANA™ Zirconia HT™ 31% - 1125 MPa		✓	✓	✓	✓	✓

*Anterior bridges only. **HT is an excellent framework material and HTML can be used for monolithic long span restorations, possibly also as a cut-back.



CEMENTATION WITHOUT CONTAMINATION

KATANA™ Cleaner removes contamination to optimise adhesive procedures. It has a high cleaning effect due to the surface active characteristic of MDP Salt. In contrast to other cleaners, KATANA™ Cleaner has a pH value of 4.5 which allows usage not only **extra-orally but also intra-orally**.



EXTRA-ORAL



INTRA-ORAL

CEMENT WITH CONFIDENCE

PANAVIA™ V5

PANAVIA™ V5 is our strongest cement ever developed. And with five shades, our most esthetic one too. This makes a totally new kind of dentistry possible, where you can be confident during cementation.



PANAVIA™ SA Cement Universal

When you require fast, easy, durable cementation to zirconia or even to glass ceramics, PANAVIA™ SA Cement Universal is the answer. Suitable for almost all indirect work, including adhesive bridges, even post placement in root canals. Easy, gum friendly excess removal, easy clean-up result, it gives you more control than ever before.



NOW, THE FINISHING TOUCH

KATANA™ ZIRCONIA TWIST DIA is suitable for polishing intra-orally delivering wonderful high-shine polishing results on zirconia thanks to its innovative shape with flexible polishing spirals.



COARSE
Pre-polishing



MEDIUM
Polishing



FINE
Highshine polishing

BORN IN JAPAN

Kuraray Europe GmbH, BU Medical Products, Philipp-Reis-Str. 4, 65795 Hattersheim am Main, Germany
+49 (0)69-30585980, dental@kuraray.com, www.kuraraynoritake.eu